TROM Advance POST-OPERATIVE BRACE





TROM Advance POST-OPERATIVE BRACE

At the push of a button, the TROM Advance represents a comprehensive post-operative solution for all patients and protocols.

Indications

- ACL, PCL, MCL & LCL surgeries
- Meniscal repairs
- Patella realignment
- Regenerative Chondroplasty
- Stable femoral fractures
- Total Knee Replacements
- High Tibial Osteotomy

Warranty: 6 months

ORDERING INFORMATION

| Part Number | Description |
|-------------|------------------------|
| 11-9114-9 | TROM Advance, Cool |
| 11-9115-9 | TROM Advance, Full |
| 11-9116-9 | TROM Advance, XL, Cool |
| 11-9117-9 | TROM Advance, XL, Full |

Recommended HCPCS Code: L1832






Ultimate Fit. Total Control.

TROM Advance leverages DonJoy's premier telescoping technology to provide the ultimate in fit, comfort and adjustability while giving clinicians a broader range of motion options.

| Feature | Benefit |
|---|---|
| Tele-Fit: Sliding mechanism allowing independent strap movement just proximal and distal of the knee | Increased strap adjustability for overall better fit Accommodates cold therapy wraps/dressings |
| Streamlined hinge (20% smaller) | Lower profileLightweight |
| Telescoping adjustability between 17 3/4" and 23 1/4" | Adapts to a variety of patients |
| Quick lock from (negative) -10° to 90° | Broad range of options for any protocol |
| Quick-release buckles | Fast and easy application |
| Cool, breathable soft goods | Improved patient comfort and aesthetics |
| Wider straps (2")Full foam version | Increased support and comfort |
| XL version | • Fits thigh circumferences from 26" to 32" |

Tele-Fit



RELATED DONJOY PRODUCTS



Iceman Cold Therapy System



Flex Knee Therapy Kit



Pain Control Device











