

# Mighty. Tough.

The improved OA Adjuster™ 3 is powered by nano MAG™ and reinforced with a tough aluminum frame making it the most durable and functional DonJoy OA knee brace.



\*\*Fleming, Beynnon, et al 2000 (Am J Sports Med. 2000;28:815-824.)





**go Day Test Drive** - Simply fit your patient in any DonJoy knee osteoarthritis brace. If they are not 100% satisfied with the brace you may return it for a full credit or refund within 90 days of brace fitting.



**EASE OF USE** - The Arthritis Foundation has honored every DonJoy knee osteoarthritis braces with its coveted "Ease-of-Use" Commendation after an independent research lab conducted rigorous testing to assure that they met the needs of people with functional limitations due to the effects of arthritis.

### **ORDERING INFORMATION PDAC ASSIGNED L1845**

Part Number	Affected Compartment	Leg
11-1590-X	Medial	Right
11-1591-X	Medial	Left
11-1592-X	Lateral	Right
11-1593-X	Lateral	Left

## Available in Graphite.

Lock and off-Load Technology™ | Floating Hinge | Swiveling D-Ring

#### **PATIENT-READY**

Brace Size	Thigh Circumference 6" Above Mid-Patella
X = 1 = XS	13 – 15.5" (33 – 39 cm)
X = 2 = S	15.5 – 18.5" (39 – 47 cm)
x = 3 = M	18.5 – 21" (47 – 53 cm)
x = 4 = L	21 – 23.5" (53 – 60 cm)
X = 5 = XL	23.5 – 26.5" (60 – 67 cm)
x = 6 = XXL	26.5 – 29.5" (67 – 75 cm)
x = 7 = XXXL	29.5 – 32" (75 – 81 cm)

#### **ACCESSORIES**

Part Number	Description	Size
25-1376-0-00000	Flexion/Extension Stop Kit	Universal
11-0075-X-00000	Cotton Lycra Suspension Sleeve	S - XXL

Call 800.336.6569 or visit www.djoglobal.com/donjoy







DJO, LLC | A DJO Global Company

**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street | Vista, CA 92081-8553 | U.S.A. DJOglobal.com/donjoy

Together in Motion.

DJO Global proudly sponsors:









