

AirLift™ PTTD Brace

Engineered to provide customized support and stabilization with simple application and adjustments, the AirLift offers the practitioner a pre-fabricated treatment modality for posterior tibial tendon dysfunction.



NEED:

- **Natural foot alignment**
- **Arch support**
- **Ankle stabilization**
- **Personalized fit**
- **Ease-of use**

A step-up in PTTD management



FEATURE:

- Adjustable arch aircell
- Cross-strap wrapped around foot and ankle
- Anatomically designed semi-rigid shells and adjustable straps
- Easy-to-use inflation Hand Bulb
- Rear entry design and simple strap adjustment

BENEFIT:

- When inflated, lifts foot arch providing secure support
- Lifts foot arch and secures brace
- Realign ankle and securely maintain desired position
- Customizes inflation of aircell for individualized support and comfort
- Promotes fast and easy application that can easily accommodate painful or inflamed feet and ankles

AirLift™ PTTD Brace

Integrated features in pre-fabricated design make the difference...

- Adjustable aircell, located under foot arch, securely lifts and holds arch in a more natural position.
- Hand Bulb custom-inflates arch aircell for personalized support and comfort, and can be used without removing foot wear for ease-of-use. (Hand Bulb included with AirLift.)
- Cross-strap wraps around foot and ankle to lift foot arch while securing brace to ankle.
- Fore foot extension securely maintains aircell position for sturdy arch support.
- Anatomically designed semi-rigid shells help realign the ankle providing secure support and stabilization.
- Adjustable straps secure brace and maintain desired ankle position for stability.
- Innovative application design allows patient to slip foot into the back of the brace, giving easy access to patients with rigid feet.



INDICATIONS:

- Posterior tibial tendon dysfunction or adult acquired flat foot
- Post-op
- Rehabilitation

ORDERING INFORMATION:

P/N		Description	Shoe Size (USA)	
Left	Right		Men	Women
02PSL	02PSR	Small AirLift PTTD Brace	Up to 7	Up to 8.5
02PML	02PMR	Medium AirLift PTTD Brace	7.5 - 11	9 - 12.5
02PLL	02PLR	Large AirLift PTTD Brace	11.5 and up	13 and up

Reference:

Can a New Brace Offer Relief for PTTD? L. Grant. *Podiatry Today*, March 2010

Effects of the AirLift PTTD Brace on Foot Kinematics in Subjects With Stage II Posterior Tibial Tendon Dysfunction Neville et al., *JOSPT* March 2009

A Radiographic Analysis Of Posterior Tibial Tendon Dysfunction Bracing, For Adult Acquired Flat Foot Deformity. Burston et al., *Presented at BOA conference*, Sept. 2008.

Solving problems presented by patients with diabetes, hindfoot varus or valgus, metatarsalgia and other conditions considered. *ORTHOPEDICS TODAY* 25: 16, 2005

Federal law restricts this device to sale by or on the order of a licensed healthcare professional.

Information subject to change.

"AirLift" is a trademark of DJO LLC, USA.

US PATENT: 6,406,450 AND OTHER PATENTS PENDING

DJO Asia-Pacific

Rm. 1905, 19/F, Tower 2, Grand Central Plaza, Shatin, Hong Kong

D: +852 310 522 37 F: +852 310 514 44 M: +852 904 825 57

crystal.chan@djoglobal.com www.DJGlobal.eu © 2010 DJO Asia-Pacific



NEVER STOP GETTING BETTER™