
# TRU-PULL LITE® PATELLOFEMORAL BRACE





# TRU-PULL LITE® PATELLOFEMORAL BRACE

### **Patient Type**

Women, men, runners, recreational athletes, and patients who want to wear a patellofemoral brace over long periods of time.

### Indication

For acute mild to moderate anterior knee pain, support for patellar instability, lateral subluxation and maltracking.

### **Features and Benefits**

 Short and powerful 8.5" length design keeps the brace at a minimal size on the leg but with maximum power to pull the patella into alignment.




Patented Dual Durometer Buttress
 System is a soft and tacky buttress
 pad which interfaces directly with the
 skin providing a consistent corrective
 force on the patella during activity.



## Powerful patella control - Done Eloquently

Tru-Pull patellofemoral braces incorporate the same patented technology of independent pull straps which apply a "true" and consistent pull on the patella during flexion and extension, thereby realigning the patella.

The new Tru-Pull Lite is DonJoy's most comfortable, low profile and shortest patellofemoral brace, yet strong enough to relieve pain!



### **ORDERING INFORMATION**

Part Number	Description	SIZE
11-0260-x	TRU-PULL LITE (Right)	XS-XXXL
11-0261-x	TRU-PULL LITE (Left)	XS-XXXL
25-9985	BUTTRESS KIT	Universal

 Sizes Available: XS–XXXL
 (Thigh circumference taken 6" above mid-patella)
 x = 5 = XL......23½" - 26½".....(60-67 cm) 

 x = 1 = XS......13" - 15½"......(33-39 cm) x = 3 = M......18½" - 21"......(47-53 cm) x = 6 = XXL.....26½" - 29½".....(67-75 cm) 

 x = 2 = S.......15½" - 18½".....(39-47 cm) x = 4 = L.......21" - 23½"......(53-60 cm) x = 7 = XXXL....29½" - 32"......(75-81 cm) 

### **ADDITIONAL TRU-PULL PRODUCTS**

- (A) TruPull Advanced System
- (B) TruPull Advanced Hinged
- (C) TruPull Wrap-a-Round
- (D) TruPull Attachment to a DonJoy Rigid Brace





1430 Decision Street, Vista, CA 92081-8553 t 800.336.6569 f 800.936.6569 www.donjoy.com © 2007 DJO, LLC












DJO provides possible coding suggestions based on publicly available information as a convenience to our customers. We do not make claims, promises or guarantees as to the availability of reimbursement for any DJO product. It is within the sole discretion of the customer to determine the appropriate billing code for a product, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payor.



# TRU-PULL LITE® PATELLOFEMORAL BRACE

### **Patient Type**

Women, men, runners, recreational athletes, and patients who want to wear a patellofemoral brace over long periods of time.

### Indication

For acute mild to moderate anterior knee pain, support for patellar instability, lateral subluxation and maltracking.

### **Features and Benefits**

 Short and powerful 8.5" length design keeps the brace at a minimal size on the leg but with maximum power to pull the patella into alignment.




 Patented Dual Durometer Buttress System is a soft and tacky buttress pad which interfaces directly with the skin providing a consistent corrective force on the patella during activity.



# Powerful patella control - Done Eloquently

Tru-Pull patellofemoral braces incorporate the same patented technology of independent pull straps which apply a "true" and consistent pull on the patella during flexion and extension, thereby realigning the patella.

The new Tru-Pull Lite is DonJoy's most comfortable, low profile and shortest patellofemoral brace, yet strong enough to relieve pain!



### ORDERING INFORMATION

Part Number	Description	SIZE
11-0260-x	TRU-PULL LITE (Right)	XS-XXXL
11-0261-x	TRU-PULL LITE (Left)	XS-XXXL
25-9985	BUTTRESS KIT	Universal

Sizes Available: XS-XXXL (Thigh circumference taken 6" above mid-patella)  $x = 5 = XL......23\frac{1}{2}$ " -  $26\frac{1}{2}$ ".....(60-67 cm) x = 1 = XS......13" - 15%2"......(33-39 cm) $x = 3 = M......18\frac{1}{2}$  - 21".....(47-53 cm) x = 6 = XXL.....261/2" - 291/2".....(67-75 cm) x = 7 = XXXL ..29½" - 32"......(75-81 cm)  $x = 2 = S \dots 15\frac{1}{2}$  -  $18\frac{1}{2}$  \dots (39-47 cm) x = 4 = L......21" - 231/2"......(53-60 cm)

## **ADDITIONAL TRU-PULL PRODUCTS**

- (A) TruPull Advanced System
- (B) TruPull Advanced Hinged
- (C) TruPull Wrap-a-Round
- (D) TruPull Attachment to a DonJoy Rigid Brace



















1430 Decision Street, Vista, CA 92081-8553 t 800.336.6569 f 800.936.6569 www.donjoy.com © 2007 DJO, LLC







