



HUBER<sup>®</sup> 360

Neuro Muscular Rehabilitation  
and Assessment



# Neuro Muscular Rehabilitation and Assessment

## EXCLUSIVE HUBER® 360 TECHNOLOGY

HUBER® 360 features the new Multiaxis Motorized Platform. It comes equipped with force sensors built into the platform and handles. Trained by a health professional, your patients will work safely on both physical and cognitive skills: that's Neuro Muscular Rehabilitation. Once mobility and flexibility have been recovered on HUBER® 360, patients can enroll into a Health Sport program in your center to maintain and improve their health.

## MULTIAXIS MOTORIZED PLATFORM™ WITH BUILT-IN FORCE SENSORS



- Integrated functional assessment adapted to all patient types
- Dynamic Posture Corrector™ for a precise and progressive workout
- Multi-directional stimulation and targeted recruitment

## SPECIFIC TRAJECTORIES FOR EACH THERAPEUTIC INDICATION



Circle Point Spiral



Seesaw Sun Dial Rosette Random



## 360° REHABILITATION

With HUBER® 360, the treatment is broken down into four movement fundamentals: Flexibility and Mobility, Dynamic Strengthening, Posture and Balance, and Resistance

### *Flexibility and Mobility*

- The Dynamic Posture Corrector™ improves stretching precision
- The Multiaxis Motorized Platform™ mobilizes joints in every plane of motion



### *Dynamic Strengthening*

- The Dynamic Posture Corrector™ insures proper positioning during exercise
- Feedback on the screen allows control of strength and direction for better recruitment

### *Posture and Balance*

- The static and dynamic force platform improves objectification of posture
- Balancing games are both fun and challenging for the patient



### *Resistance*

- Dynamic exercises adapted to the patient's physical condition increase their resistance
- Short-duration high-intensity activities for faster results

## HUBER® 360 ASSESSMENT AND PROGRESS REPORT

A real patient follow-up tool adapted to your practice.  
Guaranteed to be an easy-to-use time saver.

- Full Built in assessment
- 3 tests based on Romberg and Fukuda to measure and compare balance capacity
- 2 tests to measure area of stability and mobility restriction
- 2 tests to assess strength and coordination
- Test reports and patient progress report available
- Photo and video analysis tool
- Connect to the machine wirelessly and without a USB drive
- Send the assessment directly as a PDF to patients or prescribers
- Create your own programs and free menus
- HUBER® 360 Remote control
- Monitor patients' exercise remotely



### Size:

L=1800mm  
W= 1050mm (without guardrails)  
W=1335mm (with guardrails)  
H=2080mm

**Weight:** 274 kg (with guardrails)

### Floor size:

1,90m<sup>2</sup> (without guardrails)  
2,4 m<sup>2</sup> (with guardrails)

### Power requirements:

230V / 50-60 Hz / 2 kVA  
110V / 50-60Hz / 2kA

**Maximum platform tilt:** 10°

**Maximum platform speed:**

1 rev/s at 10° tilt

**Maximum platform load:** 150 kg

**10.4" color screen**

ref HUBER 360 MD 230V : 15-0011

ref HUBER 360 MD 110V : 15-0013

### Delivered with:

1 Samsung tablet, 1 Articulated Stool, 2 pads,  
1 repositioning block, 2 protection hand rails



articulated stool



pads



repositioning block



Samsung Tablet